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Small Volume  
Enteral Feeding Pump  
BY MOOG<sup>®</sup>

## OPERATOR'S MANUAL

**MOOG<sup>®</sup>**

Enhancing Healthcare.  
Enriching Lives.™

To Purchase, Visit [Avobus.com](http://Avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

---

Page

|                       |           |
|-----------------------|-----------|
| <b>ENGLISH</b> .....  | <b>2</b>  |
| <b>FRANCAIS</b> ..... | <b>56</b> |

---

# TABLE OF CONTENTS

Page

|   |           |
|---|-----------|
| <b>TABLE OF CONTENTS</b> .....                          | <b>2</b>  |
| <b>PRODUCT OVERVIEW</b> .....                           | <b>3</b>  |
| <b>INDICATIONS FOR USE</b> .....                        | <b>3</b>  |
| <b>LIST OF WARNINGS AND CAUTIONS</b> .....              | <b>4</b>  |
| WARNINGS .....  | 4         |
| <b>1. PUMP COMPONENTS</b> .....                         | <b>7</b>  |
| Infinity Orange Small Volume Enteral Feeding Pump ..... | 7         |
| Keypad .....  | 9         |
| Display .....   | 10        |
| Disposable Set .....                                    | 11        |
| Symbols .....   | 12        |
| <b>2. DIRECTIONS FOR USE</b> .....                      | <b>13</b> |
| Priming and Loading the Disposable Set .....            | 13        |
| Operating the Pump .....                                | 16        |
| A Single Feeding Example .....                          | 16        |
| The Battery .....                                       | 20        |
| Running on Battery Power .....                          | 20        |
| Charging the Battery .....                              | 21        |
| User Preference Settings .....                          | 22        |
| <b>3. ALARMS, MESSAGES, AND INDICATIONS</b> .....       | <b>26</b> |
| <b>4. TROUBLESHOOTING GUIDE</b> .....                   | <b>33</b> |
| <b>5. CLEANING</b> .....                                | <b>36</b> |
| <b>6. ACCESSORIES</b> .....                             | <b>37</b> |
| <b>7. SPECIFICATIONS</b> .....                          | <b>40</b> |
| <b>8. IMPORTANT PEDIATRIC CONSIDERATIONS</b> .....      | <b>42</b> |
| <b>9. ADDITIONAL TECHNICAL INFORMATION</b> .....        | <b>43</b> |
| <b>10. WARRANTY</b> .....                               | <b>47</b> |
| <b>LIST OF FIGURES</b> .....                            | <b>48</b> |
| <b>WHO TO CALL</b> .....                                | <b>51</b> |
| <b>INDEX</b> .....                                      | <b>52</b> |

To Purchase, Visit [Avobus.com](http://Avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

---

## PRODUCT OVERVIEW

The Infinity® Orange® is a rotary peristaltic enteral feeding pump designed to deliver programmed doses of enteral nutrition solutions at selectable rates. It is easy to use with a simple user interface and keypad. Moog technology allows an accuracy rate of +/- 5% and the ability to safely operate in any orientation. The Infinity Orange system includes a complete line of delivery sets that are made with materials that do not contain the plasticizer DEHP or natural rubber latex. All Infinity delivery sets feature automatic free-flow protection.

The Infinity Orange is designed for both hospital and alternate site care use. It is durable, water resistant, and easily carried in a Moog carry pack. The Infinity Orange's size, weight, accuracy, and portability promote and support health improving physical activity.

Enteral feeding is considered a medical procedure and following the instructions for use in the Operator Manual is important for the pump to operate normally.

## INDICATIONS FOR USE

The Infinity Orange enteral feeding pump is a rotary peristaltic pump designed to deliver programmed doses of enteral nutrition solutions at selectable rates.

---

## LIST OF WARNINGS AND CAUTIONS

Be sure to carefully read the following warnings and cautions. Failure to review all the warnings and cautions can cause harm to the operator or patient. Healthcare personnel who provide training to lay users should be sure to review all the warning and cautions with the lay user. If there is a change in the performance of the pump, contact customer service.

### WARNINGS

- ▶ Use ONLY feeding solutions prescribed by the responsible physician, registered dietitian, registered nurse, or other licensed practitioner. Failure to do so may result in harm to the patient.
- ▶ Use ONLY Infinity brand disposable sets to ensure proper fluid delivery. Others will not deliver the correct dose, may allow dangerous free-flow conditions, and may generate hazardous fluid pressures which may activate occlusion alarms at unpredictable pressures.
- ▶ Low flow rates, combined with high dose settings may exceed the life of the disposable set. The disposable bag and screw cap sets must be replaced every 24 hours and the spike sets must be replaced every 48 hours (max 1.5 L) to maintain delivery accuracy, allow proper air and occlusion sensing, and prevent growth of bacteria. DO NOT PROGRAM RATE AND DOSE COMBINATIONS WHICH EXCEED A 24 HOUR FEED REGIMEN FOR THE DISPOSABLE BAG AND SCREW CAP SETS OR A 48 HOUR FEED REGIMEN FOR THE SPIKE SETS.
- ▶ The battery capacity is an approximation. If you are unsure that enough capacity remains for your intended use, recharge it. Failure to do so may result in delayed or under-delivered therapy.
- ▶ To avoid electrical shock, never clean pump with charger plugged into an outlet or when pump is on.
- ▶ To avoid electrical shock, make sure the Infinity AC Adapter/Charger is completely dry before plugging into an electrical outlet.
- ▶ Do not use Infinity Orange enteral feeding pump for delivery of non-enteral solutions. Serious injury may result.
- ▶ Proper operation of pump requires door to be closed and latched. Make sure door is closed and latched when motor is running.

- 
- ▶ If an error occurs (ER01 - ER99 will appear in display), check all settings and volume counters and verify before starting a feeding. Failure to do so may result in delayed or under-delivered therapy.
  - ▶ No modification of this equipment is allowed including the pump, disposable sets, and accessories. Modification can result in delayed therapy, over or under delivery, electrocution, injury, or death.
  - ▶ Prior to starting the pump, confirm rate, dose, and food type are programmed as prescribed to avoid an unsafe feeding condition.
  - ▶ Device contains small parts. Dispose of or store connector cap when not in use to avoid choking hazard for children under 3 years.
  - ▶ Avoid leaving tubing and cords where infants, children, or those deemed at high risk for medical line entanglement, can become entangled. Be aware that if these items become wrapped around a patient's neck, there is an increased risk of strangulation or death. Caregivers of patients who are at risk for entanglement should discuss with their health care provider how to properly manage their lines and properly monitor patients based on their needs while in use.

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

---

## CAUTIONS

- ▶ Federal law (U.S.A.) restricts this device to sale by or on the order of a physician, registered dietitian, registered nurse, or other licensed practitioner.
- ▶ Dispose of the Infinity Pump, Infinity Orange Pump and Disposable Sets properly, as required by local law.
- ▶ If any leaks are detected in the disposable set, stop pump operation and disconnect set from patient and replace with a new disposable set. Failure to do so may result in under-delivered therapy.
- ▶ Care should be used when manually priming delivery set to ensure cassette is not damaged by excess force.
- ▶ Do not overload carry packs with personal items. Pump function may be affected by kinked or pinched tubing and/or unintentional button presses.
- ▶ The Infinity Orange uses an optical air-in-line sensor to detect air in the tubing as it passes through the pump, causing the pump to alarm if feeding is inadvertently interrupted or when feeding is complete. It is important to keep this sensor free of debris including spilled fluids or food that may adhere to the sensor and prevent proper detection of air in the tubing.

The clinician and caregiver should be trained regarding the care, use, and maintenance of the Infinity Orange pump and accessories, including proper cleaning of the area around the optical air-in-line sensor.

- ▶ Use ONLY commercially available pre-packed or commercially prepared feeding solutions formulated for use with a feeding pump that are prescribed by a licensed health care provider or dietitian.
- ▶ **DO NOT USE HOMEMADE BLENDERIZED OR LIQUIDIZED FOODS, OR OTHER NON-PRESCRIBED, NON-COMMERCIALLY AVAILABLE FEEDING SOLUTIONS.** Use of homemade blenderized or other non-prescribed feeding solutions may cause sticky solids which could cause an obstruction of the optical sensor. It is recommended that all users clean and maintain the pump and optical sensors according to manufacturer's recommendations.
- ▶ The use of commercially available blenderized formula (HCPCS Code: B4149) can impact pump accuracy and performance. It is recommended to follow the formula manufacture guidelines for pump use.

# 1. PUMP COMPONENTS

## Infinity Orange Small Volume Enteral Feeding Pump


Figure 1-1 Pump: Front

1. Door
2. Pump Wheel
3. Display
4. Tab
5. Upstream Pressure Sensor
6. Air Sensor
7. Downstream Pressure Sensor
8. Receptacle for Cassette
9. Keypad

## Infinity Orange Small Volume Enteral Feeding Pump

Figure 1-2 Pump: Rear

1. Serial Number Label
2. Charger/Power Connector
3. Receptacle for Pole Clamp
4. Speaker
5. Pump Instruction Label

## Keypad

 Press **ON/OFF** key for 1.5 seconds to either turn Infinity Orange on or off.

 Press **PRIME** key to rapidly fill disposable set with fluid.

Prime key is used for **hold to prime** and **auto prime**. See page 15 for hold to prime and page 23 for auto prime.

 Press **RUN/PAUSE** key to start or pause pump, or to clear an alarm.

 Press **FOOD TYPE** key to display and select type of food to be delivered (formula or human milk).

 Press **VOL/TOTAL** key once to display volume delivered for current feeding. Press key again to display total cumulative volume delivered in all feed cycles since total volume was last cleared.

 Press **RATE/DOSE** key to transition between rate and dose settings, and to confirm settings.

 Press **CLEAR** key to reset displayed function to its minimum value. Press and hold for three seconds to clear rate, dose, and volume delivered to minimum values at the same time.

 Press **+** key to increase displayed feed function. Press and hold key to increase value rapidly.

 Press **-** key to decrease displayed feed function. Press and hold key to decrease value rapidly.

NOTE: The following keys only function when pump is in pause mode: **PRIME**, **CLEAR**, **+**, and **-**.

Figure 1-3 Keypad Layout

### Display

The display includes large alphanumeric characters, as well as smaller symbols and words. All segments and symbols on the display are shown in *Figure 1-4*.

Settings and messages are displayed in large alphanumeric characters (*Figure 1-5*).

When the large alphanumeric characters display numbers, one of the smaller words below the characters indicates which function is being displayed (*Figure 1-6*).

The run symbol indicates pump is running (*Figure 1-7*).

The wall plug symbol indicates the charger is plugged in and connected to a live power outlet. The battery symbol indicates pump is running on battery power only. The fuel gauge symbol indicates how much battery charge remains (*Figure 1-8*).

NOTE: When operating on A/C power the display light will remain on. When operating on battery power the display light will automatically turn off 10 seconds after the last key is pressed. If you would prefer the light to turn off when connected to an outlet, you can change the pump light setting to OFF. See the User Preference Settings section (*page 22*) for instructions to change this setting.


Figure 1-4 Display Segments


Figure 1-5 Alphanumeric Characters and Units


Figure 1-6 Setting Indicators

Figure 1-7 Run Indicators

Figure 1-8 Power Indicators

## Disposable Set

The only disposable sets approved for use with Infinity Orange pump are:

- ▶ Infinity Orange Screw Cap Set - Order Number INFOCAP-E
- ▶ Infinity Orange 100 ml Bag Set - Order Number INF0100-E
- ▶ Infinity 500 ml Bag Set - Order Number INF0500-E
- ▶ Infinity 1200 ml Bag Set - Order Number INF1200-E
- ▶ Infinity Safety Screw Spike Set - Order Number INF0020-E

Figure 1-9 Infinity Orange  
100 ml Bag Set, or Infinity  
500 ml or 1200 ml Bag Set



Figure 1-10 Infinity  
Orange Screw Cap Set

Figure 1-11 Infinity Safety Screw Spike Set

- |  |   |
|--|---|
| 1. Cassette                                  | 6. ENFit connector                                |
| 2. "W" Symbol                                | 7. Protective cover                               |
| 3. Downstream Tubing                         | 8. Screw Cap Set                                  |
| 4. Upstream Tubing                           | 9. Protective Cover for Safety<br>Screw Spike tip |
| 5. 100 ml, 500 ml, or 1200 ml<br>Feeding Bag | 10. Infinity Safety Screw Spike                   |

To Purchase, Visit [Avobus.com](http://Avobus.com) or call 1-800-674-3655

## Symbols

Brief instructions for operation of pump as well as a brief explanation of each alarm message which pump may communicate are printed on the label attached to pump. These are not intended to be used in place of the Operator's Manual. They are simply a quick reference guide. Please read the Operator's Manual before operating the Infinity Orange.

**Consult  
Accompanying  
Documents**

*Figure 1-12 Consult  
Accompanying Documents*

The bag symbol is printed on pump door over port where upstream tubing enters pump.

*Figure 1-13 Feeding  
Bag Connection*

The patient symbol is printed on pump door over port where downstream tubing exits pump.

*Figure 1-14 Patient  
Connection*

EN 60601-1 Type BF degree of protection against electrical shock. No electrical connection to patient. Drop from any angle from height of 3 feet shall not damage pump operation.

*Figure 1-15 TYPE BF  
Shock Protection*

EN 60529 degree of protection against water entering the enclosure. Water jets from any direction shall have no effect.

**IP25**

*Figure 1-16 EN 60529*

This symbol is printed on pump delivery sets. It indicates only one patient should use each disposable set.

**SINGLE USE  
ONLY**

*Figure 1-17 Single  
Patient Use*

This symbol is printed on pump delivery sets. It indicates pump delivery sets are made with materials that do not contain the plasticizer DEHP

**This Device is  
Not Made with  
DEHP**

*Figure 1-18  
DEHP-free*

This symbol is printed on pump delivery sets. It indicates pump delivery sets should be replaced every 24 hours.

**REPLACE  
EVERY 24  
HOURS**

*Figure 1-19 Replace  
Every 24 Hours*

The Infinity pump and delivery sets are made with materials that do not contain natural rubber latex.

The Infinity pump delivery sets are made with materials that do not contain Bisphenol A (BPA) in the fluid path.

## 2. DIRECTIONS FOR USE

### Recommendation for First Use:

Since battery may not be fully charged when pump is first received, it is recommended that battery be charged for 6 hours prior to operating on battery power (see page 21 for additional information).

### Priming and Loading the Disposable Set

#### Step 1:

- ▶ **If you are using the Infinity Safety Screw Spike Set with a vented bottle or pre-filled bag:** Remove protective cover from tip, then insert tip into formula exit port of container. If using the Infinity Safety Screw Spike Set, secure the tip by tightly screwing threaded components together. (Figure 2-1).
- ▶ **If you are using the 100ml, 500ml, or 1200ml disposable set:** Hold bag upright and pour in feeding solution (Figure 2-2). Close cap securely.
- ▶ **If you are using the Screw Cap delivery set:** Fill a feeding bottle/container with a 40 mm opening and secure it to the screw cap. Place the bottle/container into the suspension bag and secure the bag to an IV pole (Figure 2-3).

NOTE: Blenderized or aggressively mixed solutions may have foam. If using this type of solution, allow it to sit for 10 to 15 minutes before pouring into bag. This will reduce the chance of an alarm due to air in the tubing.

#### Step 2:

Remove protective cover from ENFit connector (Figure 2-4).

NOTE: If set is to be used with a carry pack, all air must be removed from bag and tubing.  
**Continue to step 3 for instructions on removing air.** If set is to be hung above pump, i.e. on an IV pole, **you may skip to step 4.**

#### Step 3:


Figure 2-1 Screw Tip Into Container


Figure 2-2 Fill Bag


Figure 2-3 Fill Screw Cap Set


Figure 2-4 Remove protective cover from ENFit connector

Turn bag upside down and gently squeeze. Tilt bag as needed to keep air at tubing port (Figure 2-5).

**Step 4:**

Gently pinch orange colored tubing **below** the “” symbol. Hold this position until air is removed from tubing. Gently squeeze bag at same time to assist fluid flow. If fluid does not flow, pinch pressure may be too strong. (Figure 2-6).

NOTE: Inside the orange colored tubing, below the “” symbol is the in-line occluder. The in-line occluder is the built-in anti free-flow valve. By pinching tubing gently, the tubing moves away from the in-line occluder allowing fluid to flow (Figure 2-7). It is important to only pinch the tubing below the “” symbol to avoid damaging the in-line occluder.

NOTE: Be extremely careful to pinch tubing using only the flat, soft part of your finger. Do not use fingernails when pinching tubing; doing so may puncture the delivery set (Figure 2-8). If you suspect that the delivery set has been punctured, please discard and use a new set

NOTE: Air may also be removed from tubing using the pump priming feature. See step 7 for instructions on using the pump priming feature.


Figure 2-5 Squeeze Bag


Figure 2-6 Pinch Tubing


Figure 2-7 Tubing Segment Being Pinched


Figure 2-8 Tubing Segment Being Punctured by Fingernail

**Step 5:**

Loop silicone tubing around pump wheel stretching lightly. Place cassette into pump (Figure 2-9). **Close pump door** (Figure 2-10).

NOTE: See pages 38-39 for pump use with a backpack.

**Step 6:**

Press and hold **ON/OFF** key for 1.5 seconds to turn pump on. While pump runs through a self test, display will light and an audible alarm will sound as pump displays the nine digit serial number three digits at a time for one second each. The display will then show the letter 'N' followed by a number, which is the software revision.

Next, all segments of display will be shown for 2 seconds. **Verify all display segments and symbols are active.**

The self test is complete and pump will then display last programmed rate and will be in pause mode.

If any air is still in the tubing, continue to Step 7 to use pump prime feature.

**Step 7:**

Press and hold **PRIME** key. Alarm will sound once and pump will begin pumping at an approximate rate of 360 ml/hr. Display will read "TO STOP LET GO" (Figure 2-11). Once all air is removed from tubing, release key. Pump will stop, display will revert to last programmed rate, and pump will be in pause mode.

NOTE: occlusion and air in line alarms are disabled while priming.

NOTE: Pump also has the ability to auto prime. To enable the auto prime function, see User Preference Settings on page 23.


Figure 2-9 Seat Cassette


Figure 2-10 Close Pump Door


Figure 2-11 Priming Indication


## Operating the Pump

### A Single Feeding Example

**CAUTION: Please refer to pages vi for additional usage considerations.**

If you want to deliver 50 ml of enteral feeding solution at a rate of 100 ml/hr and then stop the pump:

Key Sequence for "A Single Feeding Example"

#### Step 1:

Prepare the disposable set with enteral feeding solution as described in the Priming and Loading the Disposable Set instructions (page 13).

#### Step 2:

Hang feeding bag or container so that the bottom of bag is at or above the level of the pump door (refer to page 42).

**OR** - If an Infinity brand carry pack is to be used, load pump and feeding bag into the proper compartments, securing pump, bag and tubing with the pack's straps (page 38).

#### Step 3:

Turn pump on by pressing the **ON/OFF** key. After running the self test, display will show last programmed rate.

Figure 2-12 Program Rate to 100 ml/hr

#### Step 4:

Press the **+** or **-** key to change the rate to 100 ml/hr. Hold down either key to change rapidly (Figure 2-12).

#### Step 5:

Press the **RATE/DOSE** key to display dose. Press the **+** or **-** key until a dose of 50 ml is displayed (Figure 2-13).


Figure 2-13 Program Dose to 50 ml

**Step 6**

Confirm food type to be delivered (formula or human milk) before starting pump by pressing the **FOOD TYPE** key. Press the **+** or **-** key to change between formula and human milk.

**Step 7**

Connect end of disposable set to the patient's enteral feeding tube (Figure 2-14).

**Step 8**

Press the **RUN/PAUSE** key. Pump will begin running. Display will show programmed rate and the arcs around the run symbol will rotate (Figure 2-15).

While pump is running, the following may be viewed by pressing the appropriate key: **RATE**, **DOSE**, and **FOOD TYPE**.

To view the amount delivered in the current feeding, press the **VOL/TOTAL** key (Figure 2-16). This counter will reset itself when the previously programmed dose has been completed and a new feeding is started, or if a feeding is interrupted and Rate or Dose is changed.

To view the amount delivered over the course of several feedings press the **VOL/TOTAL** key a second time. Display will read TOTL then the amount. This counter never resets itself. It can only be reset by the user (Figure 2-17).

While pump is running, the settings cannot be changed and the Prime feature is disabled.

If you want to stop pump at any time, press the **RUN/PAUSE** key. Display will show rate, and pump will be in pause mode. Pump will save the memory of where it was in the feeding before being paused. Or, turn pump off by pressing the **ON/OFF** key.


Figure 2-14 Connect disposable set to feeding tube

Figure 2-15 Pump is Running


Figure 2-16 Amount Currently Delivered

Figure 2-17 Total Amount Delivered

To restart pump when in paused mode, press the **RUN/PAUSE** key. The pump will restart the feeding at the point where it was paused.

NOTE: If pump is turned off and turned back on again, the volume delivered counter resets to 0 and food type resets to formula. Rate, dose, and food type must be confirmed before restarting the pump with a new feeding.

NOTE: ANY CHANGES TO PUMP SETTINGS (EXCEPT FOOD TYPE) DURING A FEEDING CYCLE WILL CAUSE PUMP TO START A NEW FEEDING; IT WILL NOT START WHERE IT LEFT OFF. RATE, DOSE, AND FOOD TYPE MUST BE CONFIRMED BY VIEWING EACH PARAMETER BEFORE RESTARTING THE PUMP.

### Feeding Completion:

When the dose has been completely delivered, pump will stop running, and display will read DOSE DONE. Pump will be in a pause mode until it is turned off or feeding is restarted. To clear DOSE DONE, press and hold the **ON/OFF** key for 1.5 seconds to turn pump off, or press the **RUN/PAUSE** key to put pump in pause mode. Adjustments to settings can be made at this time, or press the **RUN/PAUSE** key a second time to start another feeding using same rate, dose, and food type settings. Confirmation of rate, dose, and food type is required before beginning another feeding.

If you would prefer the alarm to sound when the dose is complete, you can change the pump dose complete alarm setting to BEEP WHEN DONE. See the User Preference Settings section (Page 16) for instructions to change this setting. When set to BEEP WHEN DONE, the pump will stop at the end of the feeding, beep intermittently, and display will read DOSE DONE (Figure 2-18). To silence alarm when the dose is complete, press and hold the **ON/OFF** key for 1.5 seconds to turn pump off, or press the **RUN/PAUSE** key to put pump in pause mode.

Figure 2-18 DOSE DONE

**Next Feeding:**

Turn pump on. Previous settings for rate and dose will be saved in memory. Confirm that rate, dose, and food type are correct by viewing each parameter, making necessary adjustments, and restarting the pump with a new feeding.

NOTE: If pump is turned off and turned back on again, the volume delivered counter resets to 0 and food type resets to formula. Rate, dose, and food type must be confirmed before restarting the pump with a new feeding. If dose and food type are not viewed prior to pressing the **RUN/PAUSE** key, pump will alarm and display CON FIRM DOSE AND FOOD TYPE.

NOTE: Rate, dose, and food type settings can be locked so that changes cannot be made. See the User Preference Settings section (page 22) for instructions to change these settings.

## The Battery

### Running on Battery Power

Infinity Orange will run for 24 hours at a rate of 125 ml/hr. A fully depleted battery takes approximately 6 hours to fully charge.

When pump is running on battery power the battery symbol will appear in the display (*Figure 2-19*).

*Figure 2-19 Battery Power Indicator*

The segments or blocks between E and F represent the fuel gauge of the battery. Each block represents approximately 6 hours of charge. When each bar is half spent, it will begin to flash, indicating approximately 3 hours of charge is left for that bar. When the last block is gone, the battery symbol will blink to indicate there is approximately 30 minutes of charge left. Display will flash LOW BATT every 3 seconds and pump will beep every 2 seconds to remind user of low battery condition. When battery is fully depleted pump will turn off automatically.

NOTE: The battery life is an approximation based on a fully charged battery and a flow rate of 125 ml/hr. Higher flow rates, high number of charge/discharge cycles, or extreme temperatures usage will run the battery down faster, while lower flow rates will allow the battery to last longer. The expected service life of the battery is 2 to 5 years, depending on usage.

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

## The Battery

### Charging the Battery


To charge the battery, insert plug from AC Adapter/Charger into port on the left side of pump. Plug charger into a wall outlet.

Pump will charge whether it is turned off or running.


When charger is plugged in and charging, the plug symbol will appear in the display and the fuel gauge will display four bars scrolling from E to F. This pattern continuously repeats while pump is charging (Figures 2-20a through 2-20e). When bars stop scrolling pump is done charging.

To check the status of battery during charging, disconnect charger from pump and turn on pump.


NOTE: If pump shuts down due to low battery, charging less than 10 minutes may result in an ER99 alarm. Fully charge battery if depleted (see page 27 to clear alarm).

a

b

c

d

e

Figure 2-20 Fuel Gauge

## User Preference Settings

There are **seven** user preference settings. These settings can be changed to fit the user's needs.

### Alarm Volume:

The alarm has two volumes, high or low. Pump default is low (Figure 2-21).

Figure 2-21 Alarm Volume Settings BEEP LOW and BEEP HIGH

### Settings Lock:

This allows user to set a rate, dose, and food type and then lock those settings so they cannot be changed unintentionally. This will also lock the volume counters from being unintentionally cleared. Pump default is unlocked (Figure 2-22).

Figure 2-22 Settings Lock UNLK and LOCK

### Dose Done:

There are two settings for the DOSE DONE alarm:

MUTE WHEN DONE or BEEP WHEN DONE. Pump default is MUTE WHEN DONE (Figure 2-23).

Figure 2-23 Dose Done Settings BEEP WHEN DONE and MUTE WHEN DONE

### Backlight:

To conserve battery life, the backlight will turn off automatically 10 seconds after the last key has been pressed. However, while pump is connected to an A/C current, the backlight settings can be adjusted. The light can be set to ON or OFF. With the ON setting the backlight will remain on the entire time pump is plugged in to a live power outlet. With the OFF setting, the backlight will turn on when a key is pressed and turn off 10 seconds after the last key has been pressed. Pump default is ON (Figure 2-24).

Figure 2-24 Backlight Settings LITE OFF and LITE ON

NOTE: Regardless of setting, the backlight will always turn off 10 seconds after pump is turned off.

### Max Rate:

Rate range is 0.1 ml/hr up to 300 ml/hr. However, the rate may be programmed to limit how high the rate can be set. The limits are 50, 100, 150, 200, and 250 ml/hr. Rate is not limited if user preference setting is programmed to NO MAX. Pump default is NO MAX (Figure 2-25).

Figure 2-25 MAX RATE Settings

**Max Dose:**

Dose range is 0.1 ml up to 3000 ml. However, the dose may be programmed to limit how high the dose can be set. The limits are 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, and 120 ml.

Dose is not limited if use preference setting is programmed to NO MAX. Pump default is NO MAX (Figure 2-26).

Figure 2-26 MAX DOSE Settings

**Auto Prime:**

There are two prime settings: AUTO PRIM ON or AUTO PRIM OFF (Figure 2-27). Pump default is AUTO PRIM OFF.

The auto prime feature allows the user to automatically prime the tubing set. To enable the auto prime feature, set to AUTO PRIM ON.

To operate the pump with AUTO PRIM ON:


Figure 2-27 Auto Prime Settings  
AUTO PRIM ON and  
AUTO PRIM OFF

1. Ensure tubing is disconnected from the patient.
2. Press and hold the PRIME key for 1.5 seconds.
3. Alarm will sound once and the pump will begin priming the set at an approximate rate of 360ml/hr.
4. Display will read TO STOP PUSH PRIM (Figure 2-28).
5. Once auto priming is completed, pump will stop and an alarm will sound.
6. Display will revert to the last programmed rate, and pump will be in pause mode.


Figure 2-28 TO STOP PUSH PRIM

To remove any remaining air in tubing, the user must run auto prime and stop auto prime when formula reaches the end of the set.

1. To start auto prime again, press and hold the PRIME key for 1.5 seconds.
2. Alarm will sound once, and pump will reinitiate auto prime.
3. Display will read TO STOP PUSH PRIM (Figure 2-28).
4. When the formula reaches the end of the tubing, press the PRIME button again to stop the pump.

**NOTES:**

- ▶ To stop or pause priming once the priming is initiated, press the PRIME key. If the interrupted, the next auto prime will start over, meaning it same amount of formula during each auto prime segment.


To Purchase, Visit [Avobus.com](http://Avobus.com) or call 1-800-674-3655

- ▶ When auto prime is enabled it will be the only priming feature available on the pump. To use the hold to prime only option, select AUTO PRIM OFF in the User Preference settings.

**To adjust these settings:**

Key Sequence for “Adjusting User Preference Settings”

**Step 1:**

With pump turned off, press and hold for 1.5 seconds the **+** and the **ON/OFF** keys at the **same time**. Pump will beep and then display the words BEEP HIGH or BEEP LOW. Press the **+** key to change to HIGH or the **-** key to change to LOW. Pump will give a triple beep when the **+** or **-** keys are pressed and setting will change. If setting is as desired, do not press either key.

**Step 2:**

Press the **PRIME** key. Display will read UNLK (unlocked) or LOCK (locked). Press the **+** key to change to LOCK (locked) or the **-** key to change to UNLK (unlocked). If setting is as desired, do not press either key.

**Step 3:**

Press the **PRIME** key again. Display will read MUTE WHEN DONE or BEEP WHEN DONE. Press the **+** key to change to BEEP WHEN DONE or the **-** key to change to MUTE WHEN DONE. If setting is as desired, do not press either key.

**Step 4:**

Press the **PRIME** key again. Display will read LITE ON or LITE OFF. Press the **+** key to change to LITE ON or the **-** key to change to LITE OFF. If setting is as desired, do not press either key.

**Step 5:**

Press the **PRIME** key again. Display will read MAX RATE followed by NONE or a maximum rate value such as 100 ml/hr. Press the **+** or **-** key to scroll between settings until the desired value is displayed. If setting is as desired, do not press either key.

**Step 6:**

Press the **PRIME** key again. Display will read MAX DOSE followed by NONE or a maximum dose value such as 20 ml. Press the + or - key to scroll between settings until the desired value is displayed. If setting is as desired, do not press either key.

**Step 7:** Press the **PRIME** key again. Display will read AUTO PRIM ON or AUTO PRIM OFF. Press the + key to change to AUTO PRIM ON or the - key to change to AUTO PRIM OFF. If setting is as desired, do not press either key.

**Step 8:**

Press the **ON/OFF** key to turn pump off. Setting changes will be saved automatically.

NOTE: Once user has changed the desired user preference setting, it is not necessary to view the remaining user preference settings. Each change will be automatically saved as it is made.

NOTE: Pump will turn off automatically after 2 minutes in the User Preference Settings if the user does not turn it off before then.

## 3. ALARMS, MESSAGES, AND INDICATIONS

Following is a list of all alarms, display messages, and indications that are used by the Infinity Orange enteral feeding pump. Each alarm, message, or indication is described in detail on the following pages:

| Figure             | Description  | Page |
|--------------------|--|------|
| <b>Alarms</b>      |  |      |
| 3-1                | ER01 - ER99 .....                                  | 28   |
| 3-2                | LOAD SET .....                                     | 28   |
| 3-3                | LOW BATT .....                                     | 28   |
| 3-4                | NO FLOW IN .....                                   | 29   |
| 3-5                | NO FLOW OUT .....                                  | 29   |
| 3-6                | NO FOOD .....                                      | 30   |
| 3-7                | PUSH RUN TO FEED .....                             | 30   |
| 3-8                | SHUT DOOR .....                                    | 31   |
| <b>Messages</b>    |  |      |
| 3-9                | DOSE DONE .....                                    | 31   |
| <b>Indications</b> |  |      |
| 3-10               | TO STOP LET GO .....                               | 32   |
| 3-11               | MAX .....  | 32   |
| 3-12               | CON FIRM DOSE AND FOOD TYPE .....                  | 32   |
| 3-13               | TO STOP PUSH PRIM .....                            | 33   |
|                    | Battery Doesn't Hold Its Charge .....              | 33   |
|                    | Battery Doesn't Charge .....                       | 33   |
|                    | Charger Installed but No Plug Symbol Visible ..... | 33   |

Should any of these alarms or indications continue after troubleshooting, contact your healthcare provider for pump service.

To Purchase, Visit [Avobus.com](http://Avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

## Alarms

### ER01 - ER99 (Figure 3-1)

Pump will sound a continuous beep and display will read ER followed by a two digit number to indicate a self-test has failed.

**To clear alarm:** Turn the pump off and then back on.

**WARNING:** If an error occurs (ER01 - ER99 will appear in display), all settings and volume counters should be checked and verified before starting a feeding.

### LOAD SET (Figure 3-2)

Pump will sound a dual tone beep repeatedly and display will read LOAD SET.

**Why:** This alarm will sound when pump attempted to run with the set improperly loaded or missing.

**To correct the problem:** Press the **RUN/PAUSE** key to silence alarm and place pump in pause mode, then check the following:

1. An Infinity brand disposable set is loaded into pump properly and door is closed.
2. Check cassette receptacle for cleanliness, especially around pressure sensors. If cleaning is necessary, refer to *CHAPTER 5 - CLEANING* (page 36).

### LOW BATT (Figure 3-3)

When 1 hour or less of battery life is available, display will flash LOW BATT every 3 seconds and pump will beep every 2 seconds to indicate battery is low and pump will soon stop running.

**To correct the problem:** Plug in AC Adapter/Charger to recharge battery.

Figure 3-1 ER01 - ER99 Alarm


Figure 3-2 LOAD SET Alarm


Figure 3-3 LOW BATT Alarm

### NO FLOW IN (Figure 3-4)

Pump will sound a dual tone beep repeatedly and display will read NO FLOW IN.

**Why:** This alarm has occurred because pump has detected blockage in the set between pump and bag.

**To correct the problem:** Press the **RUN/PAUSE** key to silence alarm and place pump in pause mode, then check the following:

1. Check tubing for kinks or formula clumps. Correct blockage.
2. Check cassette receptacle for cleanliness, especially around the pressure sensors. If cleaning is necessary, refer to *CHAPTER 5 - CLEANING* (page 36).

Figure 3-4 NO FLOW IN Alarm

### NO FLOW OUT (Figure 3-5)

Pump will sound a dual tone beep repeatedly and display will read NO FLOW OUT.

**Why:** This alarm has occurred because pump has detected a blockage in the set between pump and patient.

**To correct the problem:** Press the **RUN/PAUSE** key to silence alarm and place pump in pause mode, then check the following:

1. Check tubing for kinks or blockages. Remove kink or blockage.
2. Check cassette receptacle for cleanliness, especially around pressure sensors. If cleaning is necessary, refer to *CHAPTER 5 - CLEANING* (page 36).
3. Back pressure from patient may also cause this alarm. Discuss with physician.

Figure 3-5 NO FLOW OUT Alarm

**NO FOOD** (Figure 3-6)

Pump will sound a dual tone beep repeatedly and display will read NO FOOD.

**Why:** This alarm has occurred because pump has detected air in the tubing. It takes approximately 1 ml of air, which is approximately 5 inches (12.7 cm) in length entering the orange colored tubing to cause an alarm.

**To correct the problem:** Press the **RUN/PAUSE** key to silence alarm and place pump in pause mode, then check the following:

1. Check feeding bag for food. If it is empty, refill bag and prime the set to remove air from the bag and tubing.
2. Check tubing for air bubbles. If bag is full but air is present in the tubing, disconnect set from patient, press and hold the **PRIME** key until air bubble has moved to the end of the tubing.
3. Check cassette receptacle for cleanliness, especially around the air sensor. If cleaning is necessary, refer to *CHAPTER 5 - CLEANING* (page 36).
4. Check that an Infinity brand disposable set is loaded into pump properly and door is closed.
5. Check disposable set for worn tubing. If it is worn, replace with a new set.

**NOTE:** Blenderized or aggressively mixed solutions may have foam. Small foam bubbles may collect in the air sensor area and must be cleared in order to avoid a NO FOOD alarm. Allowing foamy solutions to sit for 10 to 15 minutes after mixing and prior to pouring into the bag will reduce the amount of foam.

**PUSH RUN TO FEED** (Figure 3-7)

Pump will sound a dual tone beep repeatedly and display will read PUSH RUN TO FEED.

**Why:** This alarm has occurred because pump has been in pause mode for 2 minutes.

**To correct the problem:** Press the **RUN/PAUSE** key to place pump in pause

Figure 3-6 NO FOOD Alarm

Figure 3-7 PUSH RUN TO FEED Alarm

mode again. Program and use pump or press the **ON/OFF** key to turn it off.

### **SHUT DOOR** (Figure 3-8)

Pump will sound a dual tone beep repeatedly and display will read SHUT DOOR.

**Why:** This alarm has occurred because the pump door was opened while pump was running.

**To correct the problem:** Press the **RUN/PAUSE** key to silence alarm and place pump in pause mode then shut the door. Check door is latched securely. If door will not stay closed, check door latch and latch pocket on pump are not broken or cracked.

### **Messages**

#### **DOSE DONE** (Figure 3-9)

This message indicates a feeding dose has been delivered.

**To clear pump:** Either press the **ON/OFF** key to turn pump off, or press the **RUN/PAUSE** key to put pump in pause mode before starting another feeding.

Figure 3-8 SHUT DOOR Alarm


Figure 3-9 DOSE  
DONE Message

## Indications

### TO STOP LET GO (Figure 3-10)

When pressing and holding the PRIME key during hold to prime, alarm will sound once, pump will begin pumping at an approximate rate of 360 ml/hr and display will read TO STOP LET GO.

**To stop priming:** Release the PRIME key.

### MAX (Figure 3-11)

If a rate/dose limit has been programmed in the User Preference Settings, pressing the + key when that limit has been met will cause the pump to display MAX.

For example, if the User Preference Setting for maximum rate has been programmed to 100 ml/hr when the user attempts to program a rate higher than 100 ml/hr the rate will stop increasing and each time the + key is pressed the display will read MAX. When released, display will read 100 ml/hr.

**To correct the problem:** Change limit in User Preference Settings (page 22).

### CON FIRM DOSE AND FOOD TYPE (Figure 3-12)

All feeding parameters must be confirmed before pump will begin running. This message will display if the RUN/PAUSE key is pressed prior to viewing rate, dose, and food type parameters.

**To correct the problem:** View rate, dose, and food type. Change settings where required by using the + or - key.

NOTE: The food type key allows users to select whether they are delivering formula (when selected, screen will read FORM ULA) or human milk (screen will read HUMN MILK). The pump automatically calibrates the air-in-line sensor sensitivity based on what food type setting is chosen.


Figure 3-10 Priming Indication


Figure 3-11 MAX Indication

Figure 3-12 CON FIRM DOSE AND FOOD TYPE Indication

### TO STOP PUSH PRIM (Figure 3-13)

When the auto prime is initiated, the pump screen will read TO STOP PUSH PRIM. To stop priming, push the PRIME key. See page 23 for information on auto prime.


Figure 3-13 TO STOP PUSH PRIM Indication

### Battery Doesn't Hold Its Charge

Battery will lose its charge if stored for a long period of time or if stored at high temperatures (such as in an automobile during summer).

**To correct the problem:** Recharge the battery. If the problem continues during regular use, contact your healthcare provider for service.

### Battery Doesn't Charge

If the battery symbol and the E and F of the fuel gauge (Figure 3-14) are flashing while the pump is plugged in, the battery is not charging.

Figure 3-14 Battery Symbol and the E and F of the Fuel Gauge

**To correct the problem:** Contact your healthcare provider for service.

### Charger Installed but No Plug Symbol Visible

If the plug symbol does not display after the AC Adapter/Charger is plugged into pump and a live power outlet, then charger is not charging battery.

**To correct the problem:** Check the following:

- ▶ Verify that the wall outlet works by plugging in another appliance, such as a lamp.
- ▶ Charger is connected properly.

If this does not correct the problem, pump or charger may be damaged. Contact your healthcare provider for service.

NOTE: There are no user serviceable or replaceable components inside the Infinity Orange.

## 4. TROUBLESHOOTING GUIDE

### NO FOOD Alarm

If re-priming of the pump delivery set is required, it is advisable to disconnect the enteral adapter from the patient's feeding tube while re-priming.

|   |       |   |
|---|-------|---|
| Is bag of pump delivery set empty?  | YES ▶ | Refill bag, re-prime delivery set and resume feeding.   |
| ▼ NO ▼  |       |   |
| Is air present in the tubing?   | YES ▶ | Has air been removed from <b>bag</b> of pump delivery set?  |
|   | NO ▶  | Manually prime tubing by turning bag upside down, gently squeezing lower half of bag while pinching the teal colored tubing just below the "O" drop symbol. Hold this position until all air is removed from bag and tubing.  |
|   | YES ▶ | Press and hold the <b>PRIME</b> key to move the air in the tubing past pump delivery set cassette and resume feeding.<br><br>If your pump is equipped with an auto prime option, refer to page 23 for details.  |
| ▼ NO ▼  |       |   |
| Is formula blenderized or aggressively mixed?                                     | YES ▶ | Let formula sit for 10-15 minutes before preparing pump delivery set for feeding. If foam is present in tubing, flush by priming pump delivery set either manually or with the pump and resume feeding.   |
| ▼ NO ▼  |       |   |
| Are air bubbles trapped inside the cassette?                                      | YES ▶ | Remove tubing from pump and manually prime tubing by turning bag upside down, gently squeezing lower half of bag while pinching the teal colored tubing just below the drop symbol. Hold this position until all air has passed through the teal segment of the tubing. |
| ▼ NO ▼  |       |   |
| Has the appropriate food type setting on the pump been selected for this feeding? | NO ▶  | Pause pump and change food type by pressing FOOD TYPE key and either + or - to change between formula and human milk.   |
|   | YES ▶ | If alarm continues, contact Moog Customer Service at 800.970.2337 for assistance.   |

**ER01, ER02 or ER03 Alarms**

|  |                     |  |
|--|---------------------|--|
| <b>Was pump door open while pump was turned on or while running?</b> | <b>YES/UNSURE ▶</b> | Turn pump off. Unload and reload delivery set cassette. Make sure door is closed securely, then turn pump back on. |
|  | <b>NO ▶</b>         | If alarm continues, contact Moog Customer Service at 800.970.2337 for assistance.                                  |

**NO FLOW IN or NO FLOW OUT Alarms**

- ▶ **NO FLOW IN** Occlusion between delivery set and pump.
- ▶ **NO FLOW OUT** Occlusion between pump and patient.

|  |              |   |
|--|--------------|---|
| <b>Is delivery set tubing pinched, kinked, or clogged?</b>         | <b>YES ▶</b> | Check delivery set for obstructions or kinks in tubing. Correct blockage and resume feeding.  |
| <b>▼ NO ▼</b>  |              |   |
| <b>Is tab on inside of door broken?</b>                            | <b>YES ▶</b> | Contact Moog Customer Service at 800.970.2337 to order a new door.  |
| <b>▼ NO ▼</b>  |              |   |
| <b>Is pressure sensor region of cassette receptacle clean?</b>     | <b>NO ▶</b>  | Clean pressure sensor area with a cotton swab, soft cloth, or dampened sponge; or wash entire pump under running water. Do not use abrasive materials or harsh chemicals.                 |
| <b>▼ YES ▼</b>   |              |   |
| <b>Are there visible signs of damage to pressure sensors area?</b> | <b>YES ▶</b> | Contact Moog Customer Service at 800.970.2337 to return pump for service  |
|  | <b>NO ▶</b>  | Disconnect from patient and fill a new delivery set with food. Prime set and start pump. If alarm continues, contact Moog Customer Service at 800.970.2337 to return pump for assistance. |

## LOAD SET Alarm

|  |             |  |
|--|-------------|--|
| Is door closed securely?                       | NO/UNSURE ▶ | Turn pump off. Unload and reload delivery set cassette. Make sure door is closed securely, then turn pump back on. |
| ▼ YES ▼  |             |  |
| Is door cracked or tab on side of door broken? | YES ▶       | If alarm continues, contact Moog Customer Service at 800.970.2337 to order a new door.                             |
|  | NO ▶        | If alarm continues, contact Moog Customer Service at 800.970.2337 for assistance.                                  |

## Charging Problems

|   |             |   |
|---|-------------|---|
| Is A/C adapter charger plugged into wall outlet and pump properly?          | NO/UNSURE ▶ | Check connections for both wall plug and pump adapter. Green indicator light on charger wall plug should be illuminated and pump display should have plug symbol with scrolling bars between E and F of fuel gauge. When pump is fully charged, bars will stop scrolling. |
| ▼ YES ▼   |             |   |
| Is wall outlet functioning properly?  | NO/UNSURE ▶ | Plug another device into outlet to verify outlet is functioning properly.   |
| ▼ YES ▼   |             |   |
| Does charger adapter appear to have any damage?                             | YES ▶       | Contact Moog Customer Service at 800.970.2337 to order a new charger  |
| ▼ NO ▼  |             |   |
| Does charger port appear to have any damage including bent or missing pins? | YES ▶       | Contact Moog Customer Service at 800.970.2337 to return pump for service.   |
| ▼ NO ▼  |             |   |
| Is Battery Symbol and E and F of fuel gauge flashing?                       | YES ▶       | Contact Moog Customer Service at 800.970.2337 to return pump for service.   |
|   | NO ▶        |   |

To Purchase, Visit [Avobus.com](http://Avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

## 5. CLEANING

**WARNING:** To avoid electrical shock, never clean pump or Infinity Orange AC Adapter/Charger with charger plugged into an outlet or pump turned on.

**WARNING:** Make sure the Infinity Orange AC Adapter/Charger is completely dry before plugging into an electrical outlet.

### To clean the Infinity Orange enteral feeding pump:

Pump may be cleaned with warm soapy water (standard dish soap) and a nonabrasive sponge or soft cloth. Use a cotton swab to clean pathways of cassette receptacle and to remove orange silicon residue from pump wheel roller pins.

**OR** - Apply one of the following solutions for approximately 10 minutes, then wipe pump clean with a damp cloth or sponge:

- ▶ 5% bleach and water solution
- ▶ Multipurpose household disinfectant cleaner

Rinse pump by holding under a stream of warm water. Then, dry with a clean cloth.

**NOTE:** Clean pump as needed. If the pump is used on more than one patient, clean between each patient.

**NOTE:** Avoid harsh cleaners/disinfectants. They may cause damage to pump surfaces and may affect pump functionality. Moog has tested the following disinfectants for acceptability for cleaning/disinfecting the Infinity Orange pump. These cleaning products are categorized below:

| Acceptable                            | Unacceptable  |
|---------------------------------------|---------------|
| Metrex CaviCide                       | Ball Ruthless |
| 5% Bleach and Water Solution          | Wex-Cide 128  |
| Warm Soapy Water (Standard Dish Soap) | Alcohol       |
|                                       | Vesphene Ilse |

### To clean the Carrying Packs:

The Infinity brand packs are machine-washable. Use cold water and gentle cycle, whenever possible. Hang to dry.

### To clean the AC Adapter/Charger:

Infinity Orange AC Adapter/Charger normally does not require cleaning. When desired, a dry or slightly damp cloth may be used to clean the outside surface of AC Adapter/Charger while it is disconnected from the wall outlet.


Figure 5-1 Rinsing Pump Under Stream of Water

## 6. ACCESSORIES

### Infinity AC Adapter/Charger

Order Number 23401-001

Plug AC Adapter/Charger into a wall outlet and plug connector into pump to operate on AC power and/or recharge battery (Figure 6-1).


Figure 6-1 AC Adapter/Charger

### Infinity Pole Clamp

Order Number 11981-001

Thread bolt into receptacle on the back of pump and tighten gray wheel to mount pump on clamp. If gray wheel is slightly loosened, pump can be rotated to snap into different positions. Retighten gray wheel when pump is in the desired position. Tighten black wheel to mount clamp on a pole (Figure 6-2).


Figure 6-2 Pole Clamp / Pole Clamp with Pump Attached to IV Pole

### Infinity Adjustable Angle Pole Clamp

Order Number 29152-001

Thread bolt into receptacle on the back of pump and tighten knob to mount pump on clamp. If knob is slightly loosened, pump can be rotated to snap into different positions. Retighten knob when pump is in the desired position. Tighten black wheel to mount clamp on a pole. Pull pin to adjust angle (Figure 6-3).


Figure 6-3 Adjustable Angle Pole Clamp / Adjustable Angle Pole Clamp with Pump Attached to IV Pole

### Infinity Mini Backpack

Order Numbers PCK1001 (red), PCK1002 (light green and gray), PCK1003 (black and gray)

Pack holds Infinity pump with a 100 ml or 500 ml bag in the front section or can accommodate a 1200 ml bag in the rear section. Pump is secured with a Velcro® strap. A Velcro strap secures the neck of the bag. The tubing is loaded into pump and any excess length of tubing can be secured with a Velcro tab. The downstream tubing then feeds through a port at bottom of pack (either side). Also included in this pack is a pocket that may be used to hold an ice pack.

Dimensions: 13" H x 8" W x 4.5" D. Approximate weight when loaded with pump and 100 ml delivery set fill with 100 ml of water: 2 lbs.; 500 ml delivery set filled with 500 ml of water: 2.9 lbs.; 1200 ml delivery set filled with 1200 ml of water: 4.5 lbs.

Figure 6-4 Mini Backpack


Figure 6-5 Mini Backpack (Cutaway View of Internal Compartments)


Figure 6-6 Mini Backpack (Cutaway View of Internal Straps)

### Infinity Super-Mini Backpack

Order Numbers PCK2001 (black), PCK2002 (light green and gray)

Pack holds Infinity pump with a 100 ml or 500 ml bag. Pump is secured with a Velcro strap. A Velcro strap secures the neck of the bag. The tubing is loaded into pump and any excess length of tubing can be secured with a Velcro tab. The downstream tubing then feeds through a port at bottom of pack (either side). Dimensions: 9.5" H x 8" W x 4" D. Approximate weight when loaded with pump and 100 ml delivery set fill with 100 ml of water: 1.4 lbs.; 500 ml delivery set filled with 500 ml of water: 2.3 lbs.

Figure 6-7 Super-Mini Backpack


Figure 6-8 Super-Mini Backpack (Cutaway View of Internal Compartments)


Figure 6-9 Super-Mini Backpack (Cutaway View of Internal Straps)

To Purchase, Visit [Avobus.com](http://Avobus.com) or call 1-800-674-3655

### Infinity Waist Pack

Order Numbers PCK3001 (black and gray)

Pack holds Infinity pump with a 100 ml or 500 ml bag. Pump is secured with an elastic strap. A Velcro strap secures the neck of the bag. The tubing is loaded into pump and any excess length of tubing can be secured with a Velcro tab. The downstream tubing then feeds through a port at bottom of pack (either side).

Dimensions: 7" H x 11" W x 3" D


Figure 6-10 Waist Pack


Figure 6-11 Waist Pack (Cutaway View of Internal Compartment and Straps)

### Infinity Backpack

Order Numbers PCK4001 (black and gray)

Pack holds Infinity pump with a 100 ml, 500 ml, or a 1200 ml bag. Pump is secured with a Velcro strap. A Velcro strap secures the neck of the bag. The tubing is loaded into pump and any excess length of tubing can be secured with a Velcro tab. The downstream tubing then feeds through a port at bottom of pack (either side). Also included in this pack is a pocket that may be used to hold an ice pack. Dimensions: 17" H x 8" W x 4" D.

Figure 6-12 Backpack


Figure 6-13 Backpack (Cutaway View of Internal Compartments)


Figure 6-14 Backpack (Cutaway View of Internal Straps)

## 7. SPECIFICATIONS

|                                  |  |  |
|----------------------------------|--|--|
| <b>SIZE</b>                      | Pump Dimensions:   | 1.95" H x 5.65" W x 4.05" D<br>(4.9 cm D x 14.4 cm W x 10.3 cm H)          |
|                                  | Pump Weight:   | 14.4 oz (411.0 g)  |
| <b>BATTERY</b>                   | Type:  | Lithium Ion  |
|                                  | Life:  | 24 hours @ 125 ml/hr<br>9 hours @ 300 ml/hr                                |
|                                  | Service Life:  | 2 - 5 years, depending on usage  |
|                                  | Charge Time:   | Approximately 6 hours  |
|                                  | Charge Temperature:  | 10 °C - 40 °C  |
|                                  | Additional:  | Charge Level Indicator<br>Compact Wall Charger                             |
| <b>POWER SUPPLY</b>              | AC Adapter/Charger Input:  | 100-240 VAC<br>50-60 Hz<br>400 mA  |
|                                  | Output:  | 5 VDC<br>3.0A  |
|                                  | IP Rating  | IP42   |
| <b>DISPLAY</b>                   | Backlit LCD  |  |
| <b>MODE OF OPERATION:</b>        | Continuous   |  |
| <b>OPERATING ORIENTATION</b>     | Any  |  |
| <b>FLOW RATE</b>                 | Range:   | 0.1 - 300 ml/hr  |
|                                  | Increment:   | 0.1 ml/hr from 0.1 ml/hr to 10 ml/hr<br>1 ml/hr from 10 ml/hr to 300 ml/hr |
|                                  | Accuracy:  | ± 5% or .050 ml/hr, whichever is greater<br>(see page 35)                  |
| <b>DOSE</b>                      | Range:   | 0.1 ml to 3000 ml  |
|                                  | Increment:   | 0.1 ml from 0.1 - 10 ml<br>1 ml from 10 - 3000 ml                          |
| <b>MAXIMUM INFUSION PRESSURE</b> | Maximum dead head pressure is 40 PSI   |  |
| <b>VOLUME DISPLAY</b>            | Current Dose Status<br>Accumulated   |  |
| <b>PROGRAM MEMORY RETAINED</b>   | When Pump is on or off   |  |
| <b>PEDIATRIC USE</b>             | Yes  |  |
| <b>PRIME</b>                     | Yes, hold to prime and auto prime options  |  |
| <b>DISPOSABLE SET FEATURES</b>   | Automatic Free-flow Protection<br>No Drip Chamber<br>Not made with DEHP<br>Not made with natural rubber latex<br>Expected shelf life: 3 years<br>Expiration date listed on packaging |  |

|                                   |  |   |
|-----------------------------------|--|---|
| <b>ALARMS</b>                     | ER01 - ER99  | Self-test error   |
|                                   | LOAD SET   | Set not loaded properly   |
|                                   | LOW BATT   | Battery Low   |
|                                   | NO FLOW IN   | Upstream occlusion  |
|                                   |  | -5 psi (-34 kPa) ± 3psi (21 kPa)<br>psi (pounds per square inch)<br>kPa (kilopascals)   |
|                                   | NO FLOW OUT  | Downstream occlusion,<br>12 psi (83 kPa) ± 3psi (21 kPa) and 30 sec. delay before alarm<br>15 psi (103 kPa) ± 2psi (14 kPa) and 5 sec. delay before alarm<br>18 psi (124 kPa) ± 2psi (14 kPa) and 1 sec. delay before alarm |
|                                   |  | Downstream occlusion will alarm at:<br>5 minutes at 125 ml/hr.<br>2 hours at 1 ml/hr<br>12 hours at 0.1 ml/hr   |
|                                   |  | After occlusion is cleared, up to 1 ml of fluid will be released  |
|                                   | NO FOOD  | Air in tubing - bag or set may be empty   |
| PUSH RUN TO FEED                  | Pump unused for 2 minutes  |   |
| SHUT DOOR                         | Door opened while pump was running   |   |
| <b>ACCESSORIES</b>                | AC Adapter/Charger<br>Multi-position Pole Clamp<br>Adjustable Angle Pole Clamp rotates 360°<br>Multi-purpose Infinity Custom Carry Packs |   |
| <b>WARRANTY</b>                   | 2 years  |   |
| <b>OPERATING ENVIRONMENT</b>      | Temperature:   | 41 °F to 104 °F (5 °C to 40 °C)   |
|                                   | Humidity:  | 10% to 95% non-condensing   |
|                                   | Pressure:  | 700 -1,060 hPa  |
|                                   |  | These conditions apply before and after the pump and disposable have been removed from their packaging.   |
|                                   |  | The pump requires 30 minutes to warm up/ cool down to its normal operating temperature from low/high storage temperatures.  |
| <b>STORAGE ENVIRONMENT</b>        | Temperature:   | -4 °F to 122 °F (-20 °C to 50 °C)   |
|                                   | Humidity:  | 10% to 95% non-condensing   |
|                                   |  | These conditions apply before and after the pump and disposable have been removed from their packaging.   |
| <b>SERVICE LIFE</b>               | The pump is designed to provide a minimum of 5 years of service life   |   |
| <b>CLASSIFICATION INFORMATION</b> | Infinity Orange Enteral Feeding Pump with internal battery – Type BF Applied Part<br>Infinity Orange AC Adapter / Charger – Class II     |   |

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call 1-800-674-3655

## 8. IMPORTANT PEDIATRIC CONSIDERATIONS

Infinity Orange can be used on pediatric patients if the specifications of pump meet delivery requirements of the patient. The specifications of primary importance are:

- ▶ The flow rate range of Infinity Orange is 0.1 ml/hr to 10.0 ml/hr in 0.1 ml/hr increments and 10 ml/hr to 300 ml/hr in 1 ml/hr increments.
- ▶ The pump and administration set (“pump system”) is accurate to  $\pm 5\%$  or .050 ml/hr, whichever is greater, per volume delivered under the following conditions:
  - ▶▶ utilizing Moog Infinity brand disposable sets
  - ▶▶ fluid head height at  $+6.0$  inches  $\pm 0.3$  inches ( $+15.24$  cm  $\pm 0.76$  cm) at 95% confidence and 90% reliability with respect to center of rotor.

NOTE: For rates less than 1 ml/hr, place pump level with patient.

NOTE: Head height effect on accuracy is shown below (Figure 8-1).


Figure 8-1 Head Height Effect on Accuracy

If these specifications meet the required feed regimen, Infinity Orange can be used to deliver enteral formulas to pediatric patients. **ALWAYS VERIFY RATE, DOSE, AND FOOD TYPE BEFORE PROCEEDING TO FEED.**

## 9. ADDITIONAL TECHNICAL INFORMATION

For any questions or additional information, refer to the contact information on page 44

### Infinity Orange Disposable Set Residual Fluid

When fully primed, the Infinity Orange Disposable Set tubing contains approximately 5 mL of fluid. As the Infinity Orange delivers fluid and empties the bag, approximately 3.5 mL of air is drawn into the tubing before the pump detects that the bag is empty. Be sure to account for this 1.5 mL of residual fluid when filling the disposable bag sets, especially for volume-sensitive users.

If you want to deliver 100 ml of fluid, this is how you would calculate how much fluid to place into the bag set.

| Total Fluid Desired for Feeding | + | Residual Fluid | = | Minimum Fluid Required in Bag Set |
|---------------------------------|---|----------------|---|-----------------------------------|
| 100 ml                          | + | 1.5 ml         | = | 101.5 ml                          |

Continue setting up the disposable bag set and pump according to CHAPTER 2: DIRECTIONS FOR USE, making sure to set the dose to 500 ml.

### Electromagnetic Emissions/Interference

Electromagnetic emissions may affect the operation of any electronic medical device, including enteral feeding pumps.

The Infinity Orange will not be affected by electromagnetic emissions in most environments. However, some electromagnetic fields produced by personal communication equipment, household appliances, or occupational tools may cause electromagnetic interference (EMI) which can affect the pump.

Possible sources of electromagnetic interference with electronic medical devices include, but are not limited to: cellular phones, cordless telephones, microwave ovens, anti-theft/security systems, blenders, and high-powered tools (i.e. drills, saws, chain saws). If electromagnetic emitting devices are operated within one yard/meter of the Infinity Orange, the pump may automatically shut off and settings may return to their default values. Check the pump regularly if operating near sources of electromagnetic emissions.

The Infinity Orange can safely be operated on commercial aircraft and is designed in accordance with EN 60601-1-2, EN 60601-1-4 and RTCA DO160D standards for electromagnetic emissions and immunity.


Figure 9-1 Disposable Set Residual Fluid

## 9. ADDITIONAL TECHNICAL INFORMATION

### Guidance and manufacturer's declaration – electromagnetic emissions

The Infinity Orange is intended for use in the electromagnetic environment specified below. The customer or the user of the Infinity Orange should assure that it is used in such an environment.

| Emissions Test  | Compliance | Electromagnetic environment - guidance   |
|---|------------|--|
| RF Emissions<br>CISPR 11                                | Group 1    | The Infinity Orange uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.   |
| RF Emissions<br>CISPR 11                                | Class B    |  |
| Harmonic emissions<br>IEC 61000-3-2                     | Class B    | The Infinity Orange is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/Flicker emissions<br>IEC 61000-3-3 | Complies   |  |

### Guidance and manufacturer's declaration – electromagnetic immunity

The Infinity Orange pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Infinity Orange pump should assure that it is used in such an environment.

| Immunity test  | IEC 60601 test level  | Compliance level  | Electromagnetic environment - guidance   |
|--|---|---|--|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ± 8 kV contact<br>± 15 kV air   | ± 8 kV contact<br>± 15 kV air   | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Electrostatic fast transient/burst<br>IEC 61000-4-4  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines<br>100 kHz PRR   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines<br>100 kHz PRR   | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5   | 2kV<br>(0, 90, 180, & 270° phase)   | 2kV<br>(0, 90, 180, & 270° phase)   | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <0 % $U_n$ (>95 % dip in $U_n$ ) for 0.5 cycle<br>40 % $U_n$ (60 % dip in $U_n$ ) for 5 cycles<br>70 % $U_n$ (30 % dip in $U_n$ ) for 25/30 cycles<br><0 % $U_n$ (>95 % dip in $U_n$ ) for 250/300 cycles | <0 % $U_n$ (>95 % dip in $U_n$ ) for 0.5 cycle<br>40 % $U_n$ (60 % dip in $U_n$ ) for 5 cycles<br>70 % $U_n$ (30 % dip in $U_n$ ) for 25/30 cycles<br><0 % $U_n$ (>95 % dip in $U_n$ ) for 250/300 cycles | Mains power quality should be that of a typical commercial or hospital environment. The Infinity Orange pump allows continued operation during power mains interruptions via the internal battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8   | 30 A/m  | Class 2 (30 A/m)  | Power frequency magnetic fields should be at levels characteristic of a typical location in the typical commercial or hospital environment.  |

› application of the test level.

To Purchase, Visit [Avobus.com](http://Avobus.com) or call **1-800-674-3655**

## Guidance and manufacturer's declaration – electromagnetic immunity

The Infinity Orange pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Infinity Orange pump should assure that it is used in such an environment.

| Immunity test                             | IEC 60601 test level                              | Compliance level                               | Electromagnetic environment - guidance   |
|---|---|--|--|
| Conducted RF<br>IEC 61000-4-6             | 3 Vrms<br>150 kHz to 80 MHz<br>6 Vrms (ISM bands) | 3 Vrms<br>6 Vrms (ISM bands)                   | Portable and mobile RF communications equipment should be used no closer to any part of the Infinity Orange pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br>Recommended separation distance<br>$d=1.2\sqrt{P}$  |
| Radiated RF<br>IEC 61000-4-3              | 10 V/m<br>80 MHz to 2.7 GHz                       | 10 V/m   | $d=1.2\sqrt{P}$ 80 MHz to 800 MHz<br>$d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz  |
| Radiated RF<br>RTCA/DO-160E<br>Section 20 | 75 V/m<br>100 MHz to 8 GHz                        | No Equipment<br>Category specified<br>(75 V/m) | where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Infinity Orange pump is used exceeds the applicable RF compliance level above, the Infinity Orange pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Infinity Orange pump.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile  
RF communications equipment and the Infinity Orange pump**

The Infinity Orange pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Infinity Orange pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Infinity Orange pump as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter<br>W | Separation distance according to frequency of transmitter<br>m |                                      |                                       |
|--|--|--------------------------------------|---------------------------------------|
|  | 150 kHz to 80 MHz<br>$d=1.2\sqrt{P}$                           | 80 MHz to 800 MHz<br>$d=1.2\sqrt{P}$ | 800 MHz to 2.7 GHz<br>$d=2.3\sqrt{P}$ |
| 0.01   | 0.12   | 0.12                                 | 0.23                                  |
| 0.1  | 0.38   | 0.38                                 | 0.73                                  |
| 1  | 1.2  | 1.2                                  | 2.3                                   |
| 10   | 3.8  | 3.8                                  | 7.3                                   |
| 100  | 12   | 12                                   | 23                                    |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Routine Maintenance

There is no routine calibration or adjustment procedure required for Infinity Orange. There are no user-serviceable parts. All service requests should be referred to Moog technical service.

### Disposal Information

At the end of their service lifetime, the pump and its electrical accessories (AC Adapter Charger and Data Download Accessory Cable) should be disposed of according to local standards and regulations governing the disposal of electronic waste (e-waste). Other accessories can be disposed of or recycled as standard non-hazardous waste.

### Start Up Graph and Trumpet Graph

Available upon request

## 10. WARRANTY

Solely for the benefit of the original buyer, Moog Medical Devices Group (“Moog”), warrants all new Infinity products of its manufacture (hereafter “Products”) to be free from defects in material and workmanship, and will replace or repair, F.O.B., at its factory in Salt Lake City, Utah, or other location designated by Moog, any Products returned to it within twenty-four (24) months of original purchase by the buyer. Such repair or replacement shall be free of charge.

Moog warrants to the original buyer, that the Moog-repaired portion of the Products, or replaced Products will be free from defects in material and workmanship, and Moog will replace or repair defective Products F.O.B., at its factory in Salt Lake City, Utah, or other location designated by Moog. Such Moog-performed re- pair or replacement shall carry a warranty of ninety (90) days from the date of repair or replacement or the balance of the new Product warranty as described above, whichever is greater.

THE FOREGOING WARRANTIES AND REMEDIES ARE SOLE AND EXCLUSIVE AND ARE IN LIEU OF ANY REMEDIES OR WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESS, STATUTORY OR IMPLIED. IN NO EVENT SHALL MOOG BE LIABLE FOR SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES.

This Warranty may not be modified, amended or otherwise changed, except by a written document properly executed by a duly authorized representative of Moog.

In addition, this Warranty does not apply to Products that have been altered or repaired by personnel other than those employed by Moog; nor does it apply to Products that have been subjected to misuse, abuse, neglect, improper operation of warranted Products contrary to applicable operation manuals, accident, improper maintenance or storage, Acts of God, vandalism, sabotage or fire.

MOOG PUMPS ARE DESIGNED TO BE USED ONLY WITH MOOG BRANDED ADMINISTRATION SETS. USE OF ADMINISTRATION SETS THAT HAVE NOT BEEN DESIGNED AND MANUFACTURED BY MOOG ARE CONSIDERED COUNTERFEIT PRODUCT AND WILL VOID ANY AND ALL WARRANTIES ON MOOG PUMP EQUIPMENT AND ADMINISTRATION SETS. THE USE OF COUNTERFEIT PRODUCT POSES ACCURACY AND SAFETY RISKS TO PATIENTS. CUSTOMER AGREES THAT THEY WILL NOT PURCHASE OR USE COUNTERFEIT PRODUCT WITH MOOG DEVICES.

This Warranty is void if the Product is opened or tampered with in any way without prior authorization from Moog.

This Warranty does not cover normal wear and tear and maintenance items, and specifically excludes batteries, administration sets, extension sets or any other accessory items used with the Products.

Other than Moog’s standard warranty which has been modified as stated herein, all other terms of Moog’s standard terms and conditions of sale apply, as the only terms governing this transaction.

To Purchase, Visit [Avobus.com](http://Avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

## LIST OF FIGURES

| Figure      | Description  | Page |
|-------------|--|------|
| Figure 1-1  | Pump: Front .....  | 7    |
| Figure 1-2  | Pump: Rear .....   | 8    |
| Figure 1-3  | Keypad Layout .....  | 9    |
| Figure 1-4  | Display Segments .....   | 10   |
| Figure 1-5  | Alphanumeric Characters and Units .....                                  | 10   |
| Figure 1-6  | Setting Indicators .....   | 10   |
| Figure 1-7  | Run Indicators .....   | 10   |
| Figure 1-8  | Power Indicators .....   | 10   |
| Figure 1-9  | Infinity Orange 100 ml Bag Set, or Infinity 500 ml or 1200 ml Bag Set .. | 11   |
| Figure 1-11 | Infinity Safety Screw Spike Set .....                                    | 11   |
| Figure 1-10 | Infinity Orange Screw Cap Set .....                                      | 11   |
| Figure 1-12 | Consult Accompanying Documents .....                                     | 12   |
| Figure 1-13 | Feeding Bag Connection .....   | 12   |
| Figure 1-14 | Patient Connection .....   | 12   |
| Figure 1-15 | TYPE BF Shock Protection .....   | 12   |
| Figure 1-16 | EN 60529 .....   | 12   |
| Figure 1-17 | Single Patient Use .....   | 12   |
| Figure 1-18 | DEHP-free .....  | 12   |
| Figure 1-19 | Replace Every 24 Hours .....   | 12   |
| Figure 2-1  | Screw Tip Into Container .....   | 13   |
| Figure 2-2  | Fill Bag .....   | 13   |
| Figure 2-3  | Fill Screw Cap Set .....   | 13   |
| Figure 2-4  | Remove protective cover from ENFit connector .....                       | 13   |
| Figure 2-5  | Squeeze Bag .....  | 14   |
| Figure 2-7  | Tubing Segment Being Pinched .....                                       | 14   |
| Figure 2-6  | Pinch Tubing .....   | 14   |
| Figure 2-8  | Tubing Segment Being Punctured by Fingernail .....                       | 14   |
| Figure 2-9  | Seat Cassette .....  | 15   |
| Figure 2-10 | Close Pump Door .....  | 15   |
| Figure 2-11 | Priming Indication .....   | 15   |
| Figure 2-12 | Program Rate to 100 ml/hr .....  | 16   |
| Figure 2-13 | Program Dose to 50 ml .....  | 16   |
|             | Key Sequence for "A Single Feeding Example" .....                        | 16   |
| Figure 2-14 | Connect disposable set to feeding tube .....                             | 17   |
| Figure 2-15 | Pump is Running .....  | 17   |
| Figure 2-16 | Amount Currently Delivered .....   | 17   |

|   |    |
|---|----|
| Figure 2-17 Total Amount Delivered .....  | 17 |
| Figure 2-18 DOSE DONE .....   | 18 |
| Figure 2-19 Battery Power Indicator .....   | 20 |
| Figure 2-20 Fuel Gauge .....  | 21 |
| Figure 2-23 Dose Done Settings BEEP WHEN DONE and MUTE WHEN DONE .....                                      | 22 |
| Figure 2-21 Alarm Volume Settings BEEP LOW and BEEP HIGH .....  | 22 |
| Figure 2-22 Settings Lock UNLK and LOCK .....   | 22 |
| Figure 2-24 Backlight Settings LITE OFF and LITE ON .....   | 22 |
| Figure 2-25 MAX RATE Settings .....   | 22 |
| Figure 2-26 MAX DOSE Settings .....   | 23 |
| Figure 2-27 Auto Prime Settings AUTO PRIM ON and AUTO PRIM OFF .....  | 23 |
| Figure 2-28 TO STOP PUSH PRIM .....   | 23 |
| Figure 3-1 ER01 - ER99 Alarm .....  | 27 |
| Figure 3-2 LOAD SET Alarm .....   | 27 |
| Figure 3-3 LOW BATT Alarm .....   | 27 |
| Figure 3-4 NO FLOW IN Alarm .....   | 28 |
| Figure 3-5 NO FLOW OUT Alarm .....  | 28 |
| Figure 3-6 NO FOOD Alarm .....  | 29 |
| Figure 3-7 PUSH RUN TO FEED Alarm .....   | 29 |
| Figure 3-8 SHUT DOOR Alarm .....  | 30 |
| Figure 3-9 DOSE DONE Message .....  | 30 |
| Figure 3-10 Priming Indication .....  | 31 |
| Figure 3-11 MAX Indication .....  | 31 |
| Figure 3-12 CON FIRM DOSE AND FOOD TYPE Indication .....  | 31 |
| Figure 3-13 TO STOP PUSH PRIM Indication .....  | 32 |
| Figure 3-14 Battery Symbol and the E and F of the Fuel Gauge .....  | 32 |
| Figure 5-1 Rinsing Pump Under Stream of Water .....   | 36 |
| Figure 6-1 AC Adapter/Charger .....   | 37 |
| Figure 6-2 Pole Clamp / Pole Clamp with Pump Attached to IV Pole .....                                      | 37 |
| Figure 6-3 Adjustable Angle Pole Clamp / Adjustable Angle<br>Pole Clamp with Pump Attached to IV Pole ..... | 37 |
| Figure 6-4 Mini Backpack .....  | 38 |
| Figure 6-7 Super-Mini Backpack .....  | 38 |
| Figure 6-5 Mini Backpack (Cutaway View of Internal Compartments) .....                                      | 38 |
| Figure 6-8 Super-Mini Backpack (Cutaway View of Internal Compartments) .....                                | 38 |
| Figure 6-6 Mini Backpack (Cutaway View of Internal Straps) .....  | 38 |

|  |    |
|--|----|
| Figure 6-9 Super-Mini Backpack (Cutaway View of Internal Straps) .....         | 38 |
| Figure 6-10 Waist Pack .....   | 39 |
| Figure 6-12 Backpack .....   | 39 |
| Figure 6-13 Backpack (Cutaway View of Internal Compartments) .....             | 39 |
| Figure 6-11 Waist Pack (Cutaway View of Internal Compartment and Straps) ..... | 39 |
| Figure 6-14 Backpack (Cutaway View of Internal Straps) .....                   | 39 |
| Figure 8-1 Head Height Effect on Accuracy .....                                | 42 |
| Figure 9-1 Disposable Set Residual Fluid .....                                 | 43 |

---

## WHO TO CALL

### **Additional Information:**

YOUR HEALTHCARE PROVIDER:

---

---

---

---

---

YOUR PHYSICIAN:

---

---

---

---

---

**Contact Moog Medical for additional information  
on using or maintaining the device.**

Moog Medical  
4314 ZEVEX Park Lane  
Salt Lake City, Utah 84123 USA  
Clinical and Customer Support: (800) 970-2337  
[www.infinityfeedingpump.com](http://www.infinityfeedingpump.com)

# INDEX

## Symbols

+ Key 9

- Key 9

## A

### AC Adapter/Charger 21, 37

#### Accessories

- AC Adapter/Charger 37
- Adjustable Angle Pole Clamp 37
- Backpack 39
- Mini Backpack 38
- Pole Clamp 37
- Super-Mini Backpack 38
- Waist Pack 39

#### Alarms 27

- ER01 - ER99 27
- LOAD SET 27, 83
- LOW BATT 20, 27, 83
- NO FLOW IN 28, 34, 84
- NO FLOW OUT 28, 34, 84
- NO FOOD 29, 85
- PUSH RUN TO FEED 29, 85
- SHUT DOOR 30, 86

#### Alarm Volume 22

#### Auto Prime. *See* Prime Feature

## B

### Backpack 39

#### Battery 20

- AC Adapter/Charger 21
- Battery Life 20
- Battery Symbol 20
- Charge Time 20
- Charging the Battery 21
- Fuel Gauge 20

LOW BATT Alarm 20

Plug Symbol 21

Running on Battery Power 20

## C

### Carry Packs 38

#### Cleaning

- Acceptable Cleaners/Disinfectants 36, 93
- To clean the AC Adapter/Charger 36, 93
- To clean the Carrying Packs 36
- To clean the Infinity Orange Enteral Feeding Pump 36

#### CLEAR Key 9

## D

### Directions for Use

- Removing Air 13

#### Display 9

- Battery Symbol 10
- Display Light 10
- Fuel Gauge Symbol 10
- Run Symbol 10
- Wall Plug Symbol 10

#### Disposable Set 10

- 100 ml Bag Set 11
- 1200 ml Bag Set 11
- 500 ml Bag Set 11
- Adapter Protective Cover 11, 13, 65
- Cassette 11, 65
- Downstream Tubing 11, 65
- ENFit Connector 11
- Features 40
- Infinity Safety Screw Spike tip 11
- Protective Cover for Safety Screw Spike tip 11
- Safety Screw Spike tip 11
- Upstream Tubing 11

**Disposable Set Displacement 43**  
**Dose 16, 40, 97**

## E

**Electromagnetic Emissions/  
 Interference 42, 44, 99**

**ENFit Connector 11**  
**Error, ER01 - ER99 27**

## F

**Feeding Completion 18**  
**Feeding Example (Single Feeding) 16**  
**Flow Rate 16, 42**  
**Foam 13, 29**  
**Food Type 9, 17, 18, 19, 22, 31**  
**Free-flow 40**  
**Fuel Gauge 10, 21**

## I

### Indications

Battery Doesn't Charge 32  
 Battery Doesn't Hold Its Charge 32  
 Charger Installed but No Plug  
 Symbol Visible 32  
 CON FIRM DOSE AND FOOD  
 TYPE 31, 87, 106  
 MAX 31, 87  
 TO STOP LET GO 15, 32  
 TO STOP PUSH PRIM 32, 88

### Instructions for Use

Removing Air 13

## K

### Keypad

+ Key 9

**FOOD TYPE Key 9**  
 - Key 9  
 ON/OFF Key 9  
 PRIME Key 9  
 RATE/DOSE Key 9  
 RUN/PAUSE Key 9  
 VOL/TOTAL Key 9

## M

**Maintenance 46**  
**Messages 30, 86**  
 DOSE DONE 18, 30, 86  
**Mini Backpack 38**

## O

**Occlusion (Blockage) 28**  
 ON/OFF Key 9  
**Operating Orientation 40**

## P

**Pole Clamp 8, 13, 37**  
**Prime Feature 9, 15, 23, 24, 31**  
**PRIME Key 9, 15**  
**Program Memory 17, 40**  
**Pump Components**  
 Air Sensor 7  
 Charger/Power Connector 8  
 Display 7  
 Door 7, 61  
 Downstream Pressure Sensor 7  
 Keypad 7  
 Pump Instruction Label 8  
 Pump Wheel 7, 61

Receptacle for Cassette 7  
Receptacle for Pole Clamp 8  
Serial Number Label 8  
Speaker 8  
Upstream Pressure Sensor 7

## R

Rate 16, 40  
RATE/DOSE Key 9  
Rinse 36  
RUN/PAUSE Key 9

## S

Self Test 15  
Serial Number 8, 15  
Single Feeding Example 16  
Size (Pump Specifications) 40  
Software Revision 15  
Super-Mini Backpack 38

## U

User Preference Settings 22  
Adjusting User Preference Settings 24  
Alarm Volume 22  
Auto Prime 23, 77  
Backlight 22  
Dose Done 22  
Max Dose 23  
Max Rate 22  
Settings Lock 22

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Take an interactive tour of the Infinity Pump at [infinityfeedingpump.com](http://infinityfeedingpump.com). Or, call us at **800.970.2337** to schedule a demonstration or request additional information.

Speak live with a Moog Clinical Representative for pump questions and troubleshooting guidance 24 hours per day, 7 days per week: Clinical & Customer Support **800.970.2337**

Tel: **801.264.1001**  
Toll Free: **800.970.2337**  
Fax: **801.264.1051**

4314 Zevex Park Lane  
Salt Lake City, UT 84123 USA

[infinityfeedingpump.com](http://infinityfeedingpump.com)  
[medcustomerservice@moog.com](mailto:medcustomerservice@moog.com)

For patent information, visit [www.moogmedical.com/patents](http://www.moogmedical.com/patents). Infinity® and Infinity® Orange® are registered trademarks of Zevex, Inc. Moog® is a registered trademark of Moog, Inc. ENFit® is a registered trademark of GEDSA. © 2022 Moog, Inc. Moog Medical, 4314 Zevex Park Lane, Salt Lake City, UT 84123 USA. **MADE IN THE USA** LS-28945-001 Rev. E

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