ABPM 7100

Ambulatory Blood Pressure Monitor

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To Purchase, Visit Avobus.com or call [1-800-674-3655](tel:1-800-674-3655)

901050

AMBULATORY BLOOD PRESSURE MONITOR

IEM GmbH
Gewerbepark Brand 42
52078 Aachen
Germany

Distributed by Welch Allyn

Authorized Australian Sponsor
Welch Allyn Australia (Pty) Ltd.
Unit 4.01, 2-4 Lyonpark Road
Macquarie Park, NSW 2113
Phone 1800 650 083

WelchAllyn®

Advancing Frontline Care™

EN

Initial measurement	27
24-hour measurement	28
Performing a measurement	28
Cancelling a measurement	28
Unsuccessful measurement	28
Care and Maintenance	29
Cleaning	29
Disinfection	30
Maintenance plan	31
Troubleshooting	32
Basic error sources	32
Transmission error	32
Checklist	32
Error codes	33
Limited Warranty	36
Service Policy	37
EMC Guidelines and Manufacturer Declaration	38
Patient Information - operation of the ABPM 7100	41

See www.welchallyn.com/weee

Miscellaneous symbols

	Manufacturer	 yyyy-mm	Date of manufacture
	Reference/Model number		Serial number
	Reorder/Catalog number		Batch code
	Global Trade Item Number		Protection class
	NRTL-Certification		
	Defibrillation-proof type BF applied part		MR Unsafe Poses unacceptable risks to the patient, medical staff or other persons within the MR (magnetic resonance) environment

ABPM 7100 on a step-by-step basis.
These instructions for use must be kept with the product for later use!

Clinical data

The blood pressure measuring device ABPM 7100 fulfills the requirements of the ESH (European Society of Hypertension), BHS (British Hypertension Society) and ISO 81060-2.

The device has not been tested on pregnant women, including preeclamptic patients.

EN

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

not be used for blood pressure monitoring purposes in intensive care units or during surgery!

- The ABPM 7100 must not be used in aircraft.
- The device has not been tested on pregnant women, including preeclamptic patients.

Essential Performance

The main performance features are defined as blood pressure measurement with:


- Error tolerances of the pressure gauge and measurement results are within required limits (IEC 80601-2-30).
- Maximum change value in blood pressure determination is as specified in IEC 80601-2-30.
- Cuff pressurization remains within specified limits (IEC 80601-2-30).
- An error is issued in the event that successful blood pressure measurement is impossible.

The ABPM 7100 does not issue ALARMS pursuant to IEC 60601-1-8 and is not intended for use in connection with HF surgical equipment or to clinically monitor patients in intensive care units.

Basic safety means that the patient cannot be endangered by any automatic device procedure. During any unclear conditions, the ABPM 7100 must therefore transfer into the safe Standby mode, during which the ABPM 7100 cannot automatically inflate the cuff, while this can be manually triggered by pushing the START button.

In this context, any interruption of a measurement or in automatic operation by an external influence, or the ability of the ABPM 7100 to test error conditions, is regarded as the retention or restoration of basic safety, and not as non-adherence to the main performance features.

- 3 LCD-Display
- 4 START button
- 5 DAY/NIGHT button
- 6 EVENT button
- 7 PC Interface cable port

when getting up in the morning. This individually adapts the measurement interval to the patient and assists you in the analysis of the blood pressure profile.


EVENT

The patient uses the **EVENT** button to document the time of medication or to record any events which may cause the blood pressure to rise or fall. Pressing the button will trigger a measurement, the patient should note the reason for pressing the **EVENT** button in the event log.

Warning

After an automatic measurement, you should allow at least 3 minutes to elapse before actively beginning a measurement; this will prevent longer restriction of blood circulation.

- 1 The arm cuff
- 2 Air tube
- 3 Air tube connection

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

Attention!

- Extremes of temperature, humidity or air pressure can affect measurement accuracy. Please observe the operating conditions.
- Extreme temperatures, humidity or altitude can affect the performance of the blood pressure monitor. Do not store the device near a fireplace or heating unit and do not expose it to intense sunlight. Do not place the device near a nebuliser or steam generator, as the condensation may damage it.
- The blood pressure monitor takes approx. 25 minutes to go from the minimum storage temperature of -20°C to the operating temperature of +10°C in an ambient temperature of 20°C.
- The blood pressure monitor takes approx. 25 minutes to go from the maximum storage temperature of +50°C to the operating temperature of +40°C in an ambient temperature of +20°C.

7100-10	ABPM 7100 BATT COVER REPLACEMENT	Replacement battery compartment cover
ABPM-7100CBP-UPGRA	ABPM 7100 UPGRADE KIT CBP IEM	ABPM 7100 upgrade to determine central blood pressure values (monitor serial number required)
ABPM-7100PWA-UPGRA	ABPM 7100 UPGRADE KIT PWA IEM	ABPM 7100 upgrade for pulse wave analysis (monitor serial number required)
CBP-TO-PWA-UPGRA	ABPM 7100 UPGRADE KIT CBP to PWA IEM	ABPM 7100 upgrade with CBP upgrade to PWA pulse wave analysis device (monitor serial number required)

result in impaired blood circulation.

- If the patient has limited cognitive abilities, the device may only be used under supervision.
- When used on children, the device must only be applied with special care and under permanent supervision.
- While it is still attached to a patient, the device may never be connected to a PC or other device.
- Instruct the patient to place the device in such a way that, while the cuff is inflated, the tubing is not compressed or kinked, especially during sleep.
- Petechiae, haemorrhages or subcutaneous haematoma may occur in some patients.
- Instruct the patient to turn off the device, remove the cuff, and notify the doctor if they are experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that the patient may experience some mild to moderate discomfort during a blood pressure measurement.)

45°C or under 0°C.

- Never use old, partially used non-rechargeable or rechargeable batteries together with new, unused non-rechargeable or rechargeable batteries.
- Do not attempt to recharge non-rechargeable batteries. Do not attempt to open or short-circuit non-rechargeable or rechargeable batteries. There is a risk of explosion if you do so!

Attention!

Device function

- Although zinc-carbon batteries may indicate sufficient voltage during a battery test, their output is frequently insufficient to perform 24-hour measurements. Ensure that the non-rechargeable or rechargeable batteries have sufficient power. This should be at least 2.6 V for NiMH batteries and at least 3.10 V for alkaline batteries.

	MIN	TYP	MAX
A		14,00	15,00
B		5,00	5,50
C			48,75
D	49,50		50,50

Tips on using batteries

Note

- Charge rechargeable batteries fully before first use.
- Please note that NiMH batteries only reach their full charging capacity after the 4th charge cycle.
- Recharge batteries if they have been left unused for long periods.
- Avoid deep discharge to protect your rechargeable batteries.

Attention!

Internal memory battery

- If after changing the external battery the display shows "Eboot", the internal memory battery may be empty. Please contact your dealer.

Display Segment Test	999:999 to 000:000	The display of the figures (999:999 to 000:000) is accompanied by all other symbols of the LCD in succession. Check whether all segments are correctly and fully displayed (the complete program code is checked for correctness in the background)
Current 24h time	21:45	hh:mm

If the internal test detects an error, the ABPM 7100 will indicate "E004" on the display and emit an audible signal. For safety reasons, the use of the ABPM 7100 will be locked. The faulty ABPM 7100 unit should be send back immediately for repairs to your dealer or to Welch Allyn.

EN

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

11	08:00 00:00	23:59 07:59	4 2	YES NO	YES
----	----------------	----------------	--------	-----------	-----

Setting the logs via software

To set the logs via software please refer to the respective patient management software manual.

- Note**
- Logs 1, 2 and 11 are set by default but can be changed via the patient management software.
 - Log 5 is suitable for nighttime activities (night shift).
 - Log 9 is designated as "Schellong-Test".
 - Log 11 is only available to upgraded ABPM 7100 systems in connection with HMS from version 5.0. Blood pressure measurement intervals and the 24h PWA can be set separately here. Please contact Welch Allyn for further information.

11	24 – 32 cm (9.5-12.6 in)	Adult
11L	32 – 38 cm (12.6-15.0 in)	Adult Plus
12	38 – 55 cm (15.0-21.7 in)	Large Adult

Injury.

- Examine the patient for wounds, bandages, etc.
- Question the patient regarding previous treatments.
- Observe the patient closely.
- Instruct the patient to turn off the device, remove the cuff, and notify the doctor if they are experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that the patient may experience some mild to moderate discomfort during a blood pressure measurement.)

Warning


In very rare cases materials used for and on the cuff may cause allergic reactions.

- Do not use the cuff on patients with a known hypersensitivity to epoxy resin.

Caution

Intolerances caused by the use of disinfectants.

- Wash to remove residues.
- Wash the cuff sleeve with a mild detergent in the washing machine at max. 30°C without spinning.

- Legs uncrossed
- Feet flat on the ground
- With support to the back and the arms
- With the cuff center at one level with the right atrium

Note

- During measurement, the patient should be as relaxed as possible and may not speak unless he wants to report any discomfort!
- Allow for 5 minutes of relaxation before recording the first measurement value.
- Blood pressure measurements can be influenced by the patient's position (standing, sitting, lying), by exertion or the patient's physiological state. Exclude these influencing factors to the greatest possible degree!

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

- Observe the patient closely.
- Instruct the patient to turn off the device, remove the cuff, and notify the doctor if they are experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that the patient may experience some mild to moderate discomfort during a blood pressure measurement.)

Warning

If the patient is wearing an additional ME device on the same limb for monitoring purposes, the placement and inflation of the cuff may trigger the temporary loss of the existing ME device's function.

Operation and use of the automated non-invasive blood pressure monitoring device may result in prolonged impairment to the patient's circulation or to circulation in the relevant limb.

- Examine the patient.
- Question the patient regarding previous treatments.
- Observe the patient closely.
- Instruct the patient to turn off the device, remove the cuff, and notify the doctor if they are experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that the patient may experience some mild to moderate discomfort during a blood pressure measurement.)

- The device must not be in contact with the patient during a defibrillator discharge. A discharge of this kind may damage the ABPM 7100 and cause it to display incorrect values.
The cuffs and the tube are made from non-conducting material. They therefore protect the monitor against the effects of a defibrillator discharge.
- The ABPM 7100 is not suitable for simultaneous use with HF surgical equipment.
- Measurement can be interrupted at any stage by pushing a random button. This deflates the cuff and the device can be removed.

Attention

- Do not drop the device and do not place objects on top of it.

Attention

Hygiene

- Ensure hygiene in accordance with the maintenance schedule.

may lead to motion-related artefacts or incorrect measurements. For this reason, the activity log kept by the patient must be viewed and considered in the assessment of the measurements.

- An internal reboot may occur. In this case, the device will start in the most recently used operating status.
This may be caused by internal or external influences, such as electrostatic discharges from clothing, or because the internal memory battery is empty. If the internal memory battery is empty, this error will occur when the external battery is replaced. You should therefore please contact your dealer.

Initial measurement

Note

- An initial measurement is required for starting the measurement log. The initial measurement must be checked by a physician for plausibility!

Unsuccessful measurement

1. If the display shows errors, reexamine the correct procedure during set-up and positioning of the device.
2. Dismiss the patient only after a successful manual measurement!
Inform the patient sufficiently in order to explain the situation!
3. Repeat the measurement.
4. If the display still shows errors, repeat the initial operation process.
5. For further troubleshooting measures and faults removal, please refer to the **Troubleshooting** section.

Note

- Severe malfunctions are indicated by a continuous beep.
- In the event of a continuous beep, switch off the device, remove the cuff and inform your doctor.

- It is possible to wash the cuff sleeve in the washing machine at max. 30°C. Do not spin.
- Do not use fabric softeners or other washing aids (e.g. hygiene rinses, textile deodorants). These agents can leave behind residues and damage the material.
- The cuff sleeve is not suitable for drying in a dryer.

EN

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

Except the calibration check, no further maintenance work for electronic compatibility are necessary.

have simple causes:

- Check to see that all cables are connected correctly.
- Check to see whether the ABPM 7100 and the computer are switched on.
- Check to see whether the batteries have sufficiently voltage.

Note

Some errors are combined with a continuous alarm for safety reasons. The continuous alarm can be cancelled by pressing any button. If there is residual pressure inside the cuff, open the cuff immediately.

No automatic measurements will be performed.	No manual measurements performed after application.	Valid manual measurement must always be performed after the device has been positioned.
	Incorrect log set.	Set log 1 or 2.
The measurement interval does not meet your expectations.	Incorrect log set.	The programmed log does not correspond with the set log in the ABPM 7100. Check the log manually on the device.
	No manual measurements performed after application.	Conduct manual measurement to activate the set log

+ Possible continuous alarm until a button is pressed.	Blood pressure cuff incorrectly connected.	Inspection: immediate.
	Leaky points in the cuff or the cuff tubing.	If necessary, replace the cuff.
Err 7	The memory of the blood pressure measurement device is full (a maximum of 300 measurements and events can be saved; with CBP or PWA: 260 measurements).	Delete the data in the ABP Monitor but ensure that the data was stored on your PC first.
Err 8	Measurement cancelled with a pressed button.	
Err 9	Residual pressure inside the cuff	Wait for the cuff to deflate completely.
+ Possible continuous alarm until a button is pressed.	Zero point comparison was unsuccessful.	Send the device to your specialist for inspection immediately or directly to your Welch Allyn specialist.

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

According to IEC 61000-4-11	for 1 cycle		
	70% UT for 25/30 cycles	Not applicable	
	0% UT for 250/300 cycles	Not applicable	
NOTE UT is the AC voltage before the application of the test levels.			

The ABPM 7100 is intended for operation in an electromagnetic environment as specified below. The customer or user of the ABPM 7100 should ensure it is used in such an environment.		
Immunity tests	Test level	Compliance level
Radiated disturbance variables according to IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m
Conducted disturbance variables according to IEC 61000-4-6		Not applicable

to further patient injury.

Placement and inflation of the cuff on the arm at the side of a breast amputation may lead to further injury.

- Turn off the device, remove the cuff, and notify the doctor in the event of experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that some mild to moderate discomfort may be experienced during a blood pressure measurement.)

Warning

If the patient is wearing an additional ME device on the same limb for monitoring purposes, the placement and inflation of the cuff may trigger the temporary loss of the existing ME device's function.

The operation and use of the automated non-invasive blood pressure monitoring device may lead to longer impaired blood circulation in the patient or respective limb.

- Turn off the device, remove the cuff, and notify the doctor in the event of experiencing pain, swelling, redness or numbness in the limb where the cuff is placed. (It is expected that some mild to moderate discomfort may be experienced during a blood pressure measurement.)

Attention

Damage to device

- Do not wear the ABPM 7100 while showering. If you suspect that liquid has entered the device while cleaning or using it, the device shall no longer be used on the patient.
- In the device was exposed to moisture, switch off the device and remove the batteries.
- The device may not be operated around MRI scanners or in the immediate vicinity of other medical electrical equipment.
- The device must not be in contact with the patient during a defibrillator discharge. A discharge of this kind may damage the ABPM 7100 and cause it to display incorrect values. The cuffs and the tube are made from non-conducting material. They therefore protect the monitor against the effects of a defibrillator discharge.
- The ABPM 7100 must not be used in aircraft.
- Measurement can be interrupted at any stage by pushing a random button. This deflates the cuff and the device can be removed.

Attention!

- Do not drop the device and do not place objects on top of it.

EN


To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

Correct positioning of the cuff is very important for accurate measurement, and the cuff should always be applied to the same arm.

To refit the cuff and the blood pressure monitor, please follow these instructions:

- The tube connection on the cuff must be facing upwards, see Fig.
- The cuff tube must be routed in such a way as to allow for free movement of the upper arm and should lead across your neck to the other side of your body.
- Align the cuff in such a way that no part of the cuff tube can be kinked. Align the cuff in such a way that the lower edge is about 2 cm above the bend of the elbow.
- Fit the cuff to the upper arm in such a way that you can fit one finger underneath it.
- Make sure that the artery symbol on the cuff is positioned on the arm artery (brachial artery), see Fig.
- When the cuff is correctly positioned, the metal clasp should be on the outside of the upper arm (on the elbow side). The fabric must cover the skin under the metal clasp.
- Although we recommend fitting the cuff to bare skin on the upper arm, it can also be worn over a thin shirt or blouse.
- Put on the device pouch. The length of the strap can be adjusted to allow it to be worn on the hip or on the shoulder.
- Place the ABPM 7100 in its pouch so that the cuff connection and the operating buttons are freely accessible.
- Switch on the ABPM 7100 using the **ON/OFF** button.
- Start a new blood pressure measurement by pressing the **START** button.

Warning

After an automatic measurement, you should allow at least 3 minutes to elapse before actively beginning a measurement; this will prevent longer restriction of blood circulation.

Material No.
DIR

722592
80019691 Ver. E, Revision date: 2022-02

WelchAllyn®
Advancing Frontline Care™

To Purchase, Visit [Avobus.com](https://www.avobus.com) or call [1-800-674-3655](tel:1-800-674-3655)

901050

AMBULANT BLODTRYKSAPPARAT

IEM GmbH
Gewerbepark Brand 42
52078 Aachen
Tyskland

Distribueret af Welch Allyn

Authorized Australian Sponsor
Welch Allyn Australia (Pty) Ltd.
Unit 4.01, 2-4 Lyonpark Road
Macquarie Park, NSW 2113
Phone 1800 650 083

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